

New strategy for endovascular treatment of the obstructed common femoral vein

Prof. Nils Kucher
University Hospital Zurich
Clinic for Angiology

Disclosure

Speaker name: Nils Kucher

I have the following potential conflicts of interest to report:

- Consulting: Plusmedica, Optimed, BTG
- Employment in industry
- Stockholder of a healthcare company
- Owner of a healthcare company
- Other(s): Reserach grant: Sanofi, Bayer

- I do not have any potential conflict of interest

Swiss Venous Stent Registry


Long term follow up (mean 796 days) of 121 patients treated for postthrombotic syndrome

- 35 (29%) with IVC stents
- 119 (98%) with iliac stents
- **86 (71%) stents below the inguinal ligament**
- 36 (30%) bilateral interventions


Swiss Venous Stent Registry

Stent patency at 3 years for acute iliofemoral DVT versus PTS

Primary Patency Rate




Secondary Patency Rate



Endovascular treatment of PTS


Stent thrombosis

May Thurner Anatomy



Long Rank (Mantel-Cox): $p = 0.03$

Stent below inguinal ligament



Long Rank (Mantel-Cox): $p = 0.77$



PERMANENT STENT MALFUNCTION

- 5/10 patients with stents below CVF

Issues with treatment of common femoral vein occlusions

- Often combined with postthrombotic leg inflow veins (FV, DFV)
- No dedicated CFV stent available
- Laser cut nitinol stents or woven steel alloy stents in CFV may be associated with stent malfunction, stent restenosis, stent thrombosis, stent fractures, groin discomfort
- Few and debatable results for endophlebectomy of common femoral vein

Stent fractures at and below inguinal ligament



Ideal stent for common femoral vein occlusions

- Highly flexible throughout stent length
- High radial and restitive force in areas where needed (e.g. the ligament area)
- No stent fractures
- No groing pain or discomfort
- Patency and freedom from PTS

Blueflow Venous Stent for common femoral vein (plusmedica)


First braided venous nitinol stent

- Push-Pull Delivery-System for „packing“ the stent during deployment (increase in radial force)
- Re-sheathable for allowing better deployment precision



Blueflow Venous Stent for common femoral vein (plusmedica)

- Braided meshed stent made of two 0.22mm electropolished nitinol wires
- Braiding technique with two wires enables closed loop design
- Each wire loops back when it reaches the stent tip, thus creating 14 radial force stable end loops
- Two wires are welded together at two points in the center of the stent




Blueflow case: 1

68 y, male with phlegmasia



Occlusion of the common femoral vein due to iatrogenic intramural CFV hematoma

Treatment: 12 x 100 mm Blueflow Stent



Blueflow case: 2

- 48 y, male with severe PTS and leg ulcers
- History of i.v. heroin injections CFV
- Treatment: 14 x 150 mm Blueflow Stent



Blueflow Venous Stent for common femoral vein (plusmedica)

- Two centers (Arnsberg, Zurich) initiated a post-market introduction follow-up plan in february 2018
- Overall, 36 patients were treated with the blueflow stent between february and december 2018, mainly for PTS
- No device-related complications were reported
- Overall, 19 patients had complete 3 month-follow up data in December 2018

Primary patency: 89,5% (17 of 19 patients)

Secondary patency: 94,7% (18 of 19 patients)

Summary

- Endovascular treatment of PTS often requires stent placement in the common femoral vein
- Conventional stents (laser cut nitinol or woven steel alloy) may cause problems when placed in the common femoral vein
- The blueflow stent was specifically designed for the CFV
- Initial experience shows promising results
- Further evidence is required to confirm efficacy and safety of the blueflow stent

Thank you for your attention!



New strategy for endovascular treatment of the obstructed common femoral vein

Prof. Nils Kucher
University Hospital Zurich
Clinic for Angiology